


QX200<sup>™</sup> Droplet Digital<sup>™</sup> PCR System







Droplet Digital PCR ( $ddPCR^{TM}$ ) is Bio-Rad's unique digital PCR technology. With unrivaled precision, ddPCR provides absolute quantification of target DNA or RNA molecules without the use of standard curves. ddPCR addresses the lack of scalable and practical technologies for digital PCR implementation. The new QX200 ddPCR system puts this powerful technology in your hands, ready to unveil new worlds of research at previously unattainable levels.

PCR becomes "digital" through sample partitioning and the subsequent statistical analysis of PCR product distribution across the partitions. The QX200 system uses advanced microfluidics technology to achieve partitioning on a massive scale, generating 20,000 highly uniform nanoliter-sized droplets per sample. For even higher sensitivity, combining reaction wells can produce up to millions of droplets for a single sample. Such large-scale partitioning yields a high volume of data points which, when analyzed, enable quantitative measurements at an entirely new level.

Simple and reliable, Droplet Digital PCR technology has already led to a robust set of published findings in fields such as cancer biomarker discovery, infectious diseases, genomic alterations, and gene expression. This research is only the beginning. What discoveries will you make using Droplet Digital PCR?

# AN ELEGANT AN ELEGANT EXPERIMENT PROCEDURE



Droplet Digital PCR has a simple, user-friendly experiment setup that is designed for eight samples at a time. The process easily scales up to run a 96-sample experiment in 5 hours with minimal hands-on time. When higher throughput is required, multiple 96-sample experiments can be run in a day.

Bio-Rad's QX200 ddPCR system combines water-oil emulsion droplet technology with microfluidics. The QX200 system consists of two instruments — a droplet generator and a droplet reader — and associated consumables. The droplet generator partitions each sample into 20,000 uniform nanoliter droplets containing target and background DNA in a random distribution. Each droplet serves to partition the reactions. Droplets are transferred to a 96-well PCR plate and PCR is performed to end point in a thermal cycler. The droplets stream single file through the reader for fluorescence analysis. Positive droplets, which contain at least one copy of the target DNA or RNA molecule, exhibit increased fluorescence compared to negative droplets. The fraction of PCR-positive partitions enables the target to be quantified according to the Poisson distribution.


The QX200 system works well with both TaqMan hydrolysis probes and EvaGreen fluorescence detection chemistries. In addition, its flexible design lets users optimize for either ultra-sensitive detection by combining reaction wells or higher throughput by running single wells per sample, depending on the application.




### **KEY BENEFITS**

- Achieve absolute quantification without the use of a standard curve
- Design scalable assays for high sensitivity or high throughput
- Expand applications using flexible ddPCR chemistry EvaGreen or probes









### Prepare ddPCR reaction mix

- Combine DNA sample, primers, and/or probes with one of Bio-Rad's ddPCR supermixes
- Fully validated PrimePCR ddPCR assays
   can be used

Generate droplets

- Load the ddPCR reaction mix into the wells of a droplet generator cartridge
- 8 x 20,000 droplets are generated from each run in the QX200 droplet generator
- Target DNA (—) and background DNA (—) are randomly distributed in droplets

Perform PCR with EvaGreen or hydrolysis probes

- Transfer the droplets to a 96-well PCR plate and seal the plate
- Run the PCR protocol

Read and analyze results


- After PCR, load the 96-well PCR plate into the QX200 droplet reader
- Positive and negative droplets in each sample are read
- Analyze concentrations (copies/µl) with QuantaSoft<sup>™</sup> software

### QUANTIFYING BIOMARKERS FOR CANCER



The level of sensitivity offered by the QX200 ddPCR system in quantifying cancer biomarkers overcomes the limitations posed by other methods. Using Droplet Digital PCR technology, researchers are now able to observe much finer quantitative distinctions among mutations, better pinpointing their potential roles in cancer.


Cancer-associated mutations often evade detection due to their low concentrations relative to the background of wild-type DNA in a given sample. With its high sensitivity, the QX200 system can easily scale to quantify target concentrations as low as one out of 1,000,000 (0.0001%) total copies. What was previously undetectable with other methods can now be quantified with Droplet Digital PCR.



### **SUCCESS STORY**

Detecting *T790M* mutation in epidermal growth factor receptor (*EGFR*), an important therapeutic target in some lung cancers, allows for a better understanding of resistance to tyrosine kinase inhibitor therapies. However, other techniques lacked the sensitivity to reliably identify the *T790M* allele in a high background of *EGFR* wild type. Using ddPCR, a researcher in the biotechnology industry was able to develop a highly accurate assay to quantify *T790M* concentration.

### DETECTION OF *BRAF V600E* MUTATION




High sensitivity of the QX200 system allows quantification of BRAF V600E mutation in the presence of wild-type DNA using PrimePCR ddPCR assays. A, 2-D fluorescence amplitude plot generated by QuantaSoft software shows triplicate wells of a mixed mutant:wild-type sample. The black cluster on the plot represents the negative droplets, the green cluster represents the droplets that are positive for wild-type DNA only, the blue cluster represents the droplets that are positive for mutant DNA only, and the orange cluster represents the droplets that are positive for both mutant and wild-type DNA. B, 2-D fluorescence amplitude plot shows three replicates of a wild-type-only sample. C, fractional abundance plot shows the percentage frequency (orange markers) of the mutant DNA in a wild-type DNA background. The blue markers indicate the concentration of mutant DNA (copies/µl) and the green markers indicate the concentration of wild-type DNA (copies/µl) in each of three replicate samples. All error bars generated by QuantaSoft software represent the 95% confidence interval.

## INSIGHT INTO VIRAL RESERVOIRS



Precisely quantifying viral load is crucial for characterizing disease states and developing and validating therapies. Viral reservoirs fluctuate, sometimes dropping to very low levels that are nonetheless significant. Obtaining exact and reproducible measurements of viral DNA or RNA often represents the difference between success and failure in studying many different pathogens.


The QX200 system has the sensitivity to register extremely minute quantities of viral genetic material and distinguish it from complex mixtures. Investigators can obtain actual copy numbers of viral target RNA by generating complementary DNA through reverse transcription. Along with precision and reliability, the QX200 system also offers users the throughput needed to process samples with great efficiency.



### **SUCCESS STORY**

A researcher compared quantitative PCR (qPCR) and ddPCR approaches for detecting residual human immunodeficiency virus (HIV) infection in a clinical sample set, which included an infant who was functionally cured of HIV. The researcher found that ddPCR was five times more sensitive than qPCR when measuring HIV DNA copies/million cells, and more than 20 times more accurate when measuring viral long terminal repeats (2-LTR circles).

### DETECTION OF HERPES SIMPLEX VIRUS AND β2 MICROGLOBULIN




ddPCR enables precise and reproducible detection of herpes simplex virus 1 (HSV-1) and ß2 microglobulin (B2M) targets in a duplex assay. A, 2-D fluorescence amplitude plot shows three merged replicate samples of 10,000 copies/well of HSV-1 duplexed with 20,000 copies/well of B2M. The black cluster on the plot represents the negative droplets, the green cluster represents the droplets that are positive for HSV-1 only, the blue cluster represents the droplets that are positive for B2M only, and the orange cluster represents the droplets that are positive for both HSV-1 and B2M targets. B, concentration plot shows merged triplicate wells across an HSV-1 sample dilution series from 0 to 40,000 copies in a constant background of 20,000 copies (66 ng human genomic DNA) of B2M. The blue markers indicate HSV-1 copies/µl and the green markers indicate B2M copies/µl. All error bars generated by QuantaSoft software represent the 95% confidence interval.

### DISTINGUISHING GENOMIC VARIATIONS



Copy number variation (CNV) is a prominent source of interindividual variability in the human genome. CNV has also been associated with cancers, neurological and autoimmune diseases, and adverse drug responses. Reliably identifying such CNVs represents an important capability in cutting-edge research.


The major technical hurdle in copy number (CN) assessment is the ability to discriminate, with statistical confidence, between consecutive CN states. Fundamentally, as CN state increases, the percentage difference in target genomic material between states decreases, making it harder to measure CNV at higher orders. By partitioning the CNV amplification reaction across thousands of droplets, the QX200 system can resolve high consecutive copy number states (such as five vs. six) with an excellent degree of statistical confidence.



### **SUCCESS STORY**

A researcher investigating a highly duplicated gene involved in human brain size evolution needed a method to validate small CN variations in test samples. Droplet Digital PCR proved much more effective than other methods the researcher had previously used, showing greater precision in detecting small differences at high CN states, as well as excellent reproducibility. These improvements in CNV quantification will allow the researcher to ask more complex questions regarding neurogenetic disorders such as micro/macrocephaly, schizophrenia, and autism.

### DETECTION OF COPY NUMBER FOR MYELOCYTOMATOSIS ONCOGENE AND $\beta$ -GLOBIN



The QX200 system provides superior resolution with copy number (CN) calls from CN1 to CN13. A, 2-D fluorescence amplitude plot shows duplicate wells of a CN6 sample duplexed with myelocytomatosis oncogene (*MYC*) and β-globin (*HBB*) assays. The black cluster on the plot represents the negative droplets, the green cluster represents the droplets that are positive for *MYC* only, the blue cluster represents the droplets that are positive for *HBB* only, and the orange cluster represents the droplets that are positive for both *MYC* and *HBB* targets. B, copy number plot shows samples ranging from CN1 to CN13 with precise duplicate well values. All error bars generated by QuantaSoft software represent the 95% confidence interval. NTC, no template control.

## VALIDATION AND LIBRARY QUANTIFICATION



Bio-Rad's ddPCR platform can be easily incorporated into the next-generation sequencing (NGS) library preparation workflow to precisely quantify and balance sequencing libraries on sequencers. In addition to precise quantification, the data plots generated by the QX200 system are rich with qualitative library information, which is not available with other currently available methodologies. Digital PCR enables consistent library loading and efficient utilization of next-generation sequencers. After the NGS run, ddPCR can be used to validate the results obtained by NGS, including genomic alterations such as single nucleotide polymorphisms, mutations, and copy number variations.

### GENE EXPRESSION STUDIES

With its high sensitivity, the QX200 system is especially well-suited for discerning and quantifying rare transcripts, offering a better understanding of actual RNA function. With its high precision, the QX200 system can be used to detect small fold changes in expression of a target gene between samples, providing insights into the gene's physiological implications. The ability to measure absolute values of RNA molecules adds a new level of insight to gene expression studies. The QX200 ddPCR system allows for such measurements without requiring a standard curve.



### **SUCCESS STORY**


A research group compared ddPCR with other approaches as a method for checking Illumina TruSeq library concentrations. Using ddPCR quantification with five libraries, each with a different concentration, the researchers were able to achieve the optimal cluster density of approximately 800,000/mm<sup>2</sup> for performing NGS on each library using the MiSeq platform. Tests for library balancing showed that libraries could be balanced to within less than 15% of each other with a confidence interval of 95%.



### **SUCCESS STORY**

A researcher used ddPCR to study the differential expression of huntingtin (HTT) mRNA encoded by both wild-type and mutant HTT genes. Quantifying these different transcripts of RNA gave the researcher more insight into their cellular functions, and laid the groundwork for better understanding disease mechanisms.

### VALIDATION FOR NEXT-GENERATION SEQUENCING




02.0

Gene	ddPCR, copies/well		MiSeq, RPKM			
	100 ng RNA 4 replicates	1,000 ng 4 replicates	100 ng 2 replicates		1,000 ng 2 replicates	
GAPDH	1,671 ± 115*	16,275 ± 479*	974	953	1,077	1,061
В2М	504 ± 46*	3,450 ± 155*	233	251	229	237
SDHA	139 ± 14*	1,131 ± 81*	51	50	59	65
HPRT1	15,781 ± 2,310	140,705 ± 11,059	34	22	25	25
TBP	3,650 ± 178	31,625 ± 1,010	3	3	6	3
GUSB	1,794 ± 53	15,731 ± 1,134		12	13	12

ddPCR provides nearly 1,000-fold more sensitivity and an amplification bias–free approach compared to RNA-Seq for validation. A, concentration plot shows a twofold dilution of cDNA generated using Bio-Rad's iScript™ advanced cDNA synthesis kit for RT-qPCR. Two independent measurement sets were made: one at a high concentration range and one at a lower range, with four points overlapping. B, table shows comparative results between ddPCR and MiSeq. ddPCR detected thousands of copies per well with low-abundance transcripts (TBP and GUSB genes) using an input of 100 ng total RNA, which demonstrates enhanced sensitivity. The MiSeq sequencer detected only single-digit reads per kilobase per million reads (RPKM) for these transcripts. ddPCR, Droplet Digital PCR.

### GENE EXPRESSION ANALYSIS OF STEM CELL MARKERS



High sensitivity and precision of the QX200 system allow reliable detection of small fold changes with rare transcripts. The  $TGF\beta$  signaling pathway plays a role in ovarian stem cell regulation, and the addition of  $TGF\beta$  will increase the number of stem cells in the population, increase the expression of putative stem cell markers, and decrease terminal differentiation markers like E-cadherin (ECAD). Data show expression levels of cKIT, CD133, OCT3/OCT4, and ECAD genes across ten samples with and without  $TGF\beta$  treatment, normalized to PPIA in this study. Reproducible detection of cDNA copies as low as 0.5–1 target/µl in a 20 µl reaction was achieved. S, sample.

Bio-Rad Laboratories, Inc.

### ABSOLUTE QUANTIFICATION WITH EVAGREEN



Through droplet partitioning and statistical analysis, Droplet Digital PCR provides absolute quantification without the need for a standard curve, bringing substantial improvements to many applications that require detection of low-abundance target copies or distinction of small fold changes between samples. Combining ddPCR with the power of EvaGreen detection chemistry enables genotyping and gene expression analysis for low-level target genes. The QX200 ddPCR system enables absolute quantification using EvaGreen chemistry, making such assays ideal for target DNA measurement, viral load analysis, and microbial quantification.

### VALIDATED ddPCR ASSAYS

Bio-Rad's PrimePCR ddPCR assays are fully wet-lab validated assays, expertly designed to provide single-copy PCR resolution without the use of a standard curve. Researchers can combine the power of ddPCR with the guaranteed performance of these assays for applications involving mutation detection and copy number variation. These are the only assays available that have been validated for digital PCR.

The ease of use and reliable results these ddPCR assays provide will accelerate discovery and optimize research strategies.



### **SUCCESS STORY**


Precision cancer medicine involves rapid assessment of cancer mutations, with certain genetic errors serving as indicators for a specific therapy. A researcher was challenged with the composition and degradation of clinical samples that were used for assessment of FLT3 mutations for acute myelogenous leukemia. By partitioning the samples' DNA templates, ddPCR removed the effects of PCR inhibitors and bias in the data, resulting in a more reliable quantification of the mutations present in the samples.



### **SUCCESS STORY**


Designing an optimal probe assay may pose challenges and delay generation of reliable data. A researcher spent a great deal of time trying to design probes and primers for *EGFR L858R*, a receptor for epidermal growth factor, which is associated with lung tumors. Using the fully validated PrimePCR ddPCR mutation detection probe assay for *EGFR L858R* eliminated these tedious optimization steps, allowing the researcher to focus on the biology of the investigation, rather than assay design.

### ABSOLUTE QUANTIFICATION USING EVAGREEN



ddPCR enables precise and reproducible detection calls across a range of sample concentrations. A, concentration plot shows data for both merged wells (solid blue markers) and individual replicate wells (open blue markers) across a dilution series. Interleukin 1 (*IL-1*) samples ranging from 0 to 50,000 copies were detected with an EvaGreen assay and plotted on a log scale. There is excellent linearity and reproducibility between replicates. B, 1-D fluorescence amplitude plot shows a sample dilution curve detected with an EvaGreen *IL-1* assay in the FAM/EvaGreen channel of the QX200 droplet reader. The data demonstrate good separation between positive ( ) and negative ( ) data points. All error bars generated by QuantaSoft software represent the 95% confidence interval. NTC, no template control.

### DETECTION OF EGFR L858R MUTATION



PrimePCR ddPCR assays are fully validated for the detection of EGFR L858 mutation. A, 2-D fluorescence amplitude plot shows triplicate wells of a mixed mutant:wild-type sample. The black cluster on the plot represents the negative droplets, the green cluster represents the droplets that are positive for wild-type DNA only, the blue cluster represents the droplets that are positive for mutant DNA only, and the orange cluster represents the droplets that are positive for both mutant and wild-type DNA. B, fractional abundance plot shows the percentage frequency (orange markers) of the mutant DNA in a wild-type DNA background. The blue markers indicate the concentration of mutant DNA (copies/µl) and the green markers indicate the concentration of wild-type DNA (copies/µl) in each of three replicate samples. All error bars generated by QuantaSoft software represent the 95% confidence interval.

Bio-Rad Laboratories, Inc.



### **More Droplet Digital PCR Resources**

Specifications			Catalog #	Description			
QX200 Drople	et Generator		186-3023	Droplet PCR Supermix, 2 ml (2 x 1 ml), 200 x 20 µl reactions,			
Starting sample size 20 µl			2x supermix				
Capacity		1-8 samples/cartridge	186-3024	Droplet PCR Supermix, 5 ml (5 x 1 ml), 500 x 20 µl reactions,			
· · ·		20,000		2x supermix			
Dimensions (W x D x H) 28 x 36 x 13 cm (11 x 14 x 5")		186-3025	Droplet PCR Supermix, 25 ml (5 x 5 ml), 2,500 x 20 µl reactions,				
QX200 Droplet Reader				2x supermix			
Precision ±10%			186-3021	One-Step RT-ddPCR Kit for Probes, 2 ml (2 x 1 ml),			
		5 orders of magnitude		200 x 20 µl reactions, 2x RT-ddPCR mix, includes			
, ,		1–96 samples		1 manganese acetate tube			
Droplets per 96-well plate		Approximately 1,500,000	186-3022	One-Step RT-ddPCR Kit for Probes, 5 ml (5 x 1 ml),			
Sample illumination		Light-emitting diodes		500 x 20 μl reactions, 2x RT-ddPCR mix, includes			
Sample detection		Multipixel photon counter		2 manganese acetate tubes			
Detection char		FAM (EvaGreen), HEX (VIC)	186-4033	QX200 ddPCR EvaGreen Supermix, 2 ml (2 x 1 ml),			
Dimensions (W		66 x 52 x 29 cm (26 x 20 x 11")		200 x 20 µl reactions, 2x supermix			
	,	,	186-4034	QX200 ddPCR EvaGreen Supermix, 5 ml (5 x 1 ml),			
Ordering Information			100 1005	500 x 20 µl reactions, 2x supermix			
_			186-4035	QX200 ddPCR EvaGreen Supermix, 25 ml (5 x 5 ml),			
Catalog # Description		100 1000	2,500 x 20 µl reactions, 2x supermix				
QX200 Droplet Digital PCR System			186-4036	QX200 ddPCR EvaGreen Supermix, 50 ml (10 x 5 ml),			
186-4001	QX200 Droplet Digital PCR System, includes droplet generator,		186-3052	5,000 x 20 µl reactions, 2x supermix ddPCR Buffer Control Kit for Probes, 9 ml (2 x 4.5 ml), 2x buffer			
	droplet reader, laptop computer, software, associated component		186-4052	ddPCR Buffer Control Kit for EvaGreen, 9 ml (2 x 4.5 ml), 2x buffer			
consumables			Thermal Cycler and Plate Sealer				
186-4002	QX200 Droplet Generator, includes droplet generator,		185-1196	C1000 Touch™ Thermal Cycler with 96-Well Fast Reaction			
	1 pkg of 24 DG8 o	cartridges, 1 pkg of 24 DG8 gaskets,	103-1190	Module, includes C1000 Touch thermal cycler chassis, 96-well fast			
	2 cartridge holder	s, 1 power cord		reaction module, USB flash drive			
186-4003	•	Reader, includes droplet reader, ddPCR manual,	181-4000	PX1 <sup>™</sup> PCR Plate Sealer, includes heat sealing instrument, plate			
	2 plate holders, U	SB cable, power cord	101 4000	support block that holds 96-well and 384-well plates, sealing frame,			
186-4007				power cord			
		ges, 5 pkg of 24 DG8 gaskets		power cord			
186-4008	_	for QX100 <sup>™</sup> /QX200 Droplet Generator,	F O !	tundens du et Dietime Inn Die Deutschenden Inn in lieuwend			
	1 pkg of 24 cartrid	9	EvaGreen is a trademark of Biotium, Inc. Bio-Rad Laboratories, Inc. is licensed				
186-3009		QX100/QX200 Droplet Generator,		c. to sell reagents containing EvaGreen dye for use in real-time PCR, for			
100 0051	1 pkg of 24 gaske		research purposes only. FAM and VIC are trademarks of Applera Corporation. Illumina,				
186-3051	DG8 Cartridge H		MiSeq, and TruSeq are trademarks of Illumina, Inc. TaqMan is a trademark of Roche Molecular Systems, Inc.				
186-3030	•	or Oil for Probes, 2 x 7 ml					
186-3005 186-4005	•	or Oil for Probes, 10 x 7 ml	This product and/or its use is covered by claims of U.S. patents, and/or pending				
186-4006	•	or Oil for EvaGreen, 2 x 7 ml or Oil for EvaGreen, 10 x 7 ml	U.S. and non-U.S. patent applications owned by or under license to Bio-Rad				
186-3031	Droplet Reader		Laboratories, Inc. Purchase of the product includes a limited, non-transferable right				
186-3004	•	et Reader Oil, 1 X 1 L		under such intellectual property for use of the product for internal research purposes			
ddPCR Reag	•	5H, Z X I L	, ,	s are granted for diagnostic uses. No rights are granted for use of			
186-3026 ddPCR Supermix for Probes, 2 ml (2 x 1 ml), 200 x 20 µl reactions,				or commercial applications of any kind, including but not limited to g, quality control, or commercial services, such as contract services or			
100 0020	2x supermix		,	es. Information concerning a license for such uses can be obtained from			
186-3010	ddPCR Supermix for Probes, 5 ml (5 x 1 ml), 500 x 20 µl reactions,			pratories. It is the responsibility of the purchaser/end user to acquire any			
.00 0010	2x supermix			illectual property rights that may be required.			
186-3027	'	ddPCR Supermix for Probes, 25 ml (5 x 5 ml),					
	•	2,500 x 20 µl reactions, 2x supermix		I-time thermal cyclers are covered by one or more of the following			
186-3028	ddPCR Supermix for Probes, 50 ml (10 x 5 ml),			or their foreign counterparts owned by Eppendorf AG: U.S. Patent			
	5,000 x 20 µl reactions, 2x supermix		Numbers 6,767,512 and 7,074,367.				



### Москва

ул. Космонавта Волкова, 10, Тел./факс (495) 745-0508 sales@dia-m.ru

### Казань

Оренбургский тракт, 20, оф. 217 пер. Семашко, 114 Тел./факс (843) 277-6040 kazan@dia-m.ru

### Новосибирск

пр. Ак. Лаврентьева, 6/1 Тел./факс (383) 328-0048 nsk@dia-m.ru

### Ростов-на-Дону

Тел./факс (863) 250-0006 rnd@dia-m.ru

### Санкт-Петербург

ул. Профессора Попова, 23 Тел./факс (812) 372-6040 spb@dia-m.ru

### Пермь

Представитель в УФО -Арсения Шелемба Тел./факс (342) 202-2239