

Techfors

The compact bioreactor for laboratory, pilot plant and production







A new emphasis on instrumentation and construction

Innovative

Techfors uses new approaches to old problems, e.g. a fluid-cooled motor with no need for a circulation pump for a very quiet running of the machine.

Adaptable

Techfors covers a range of vessel sizes up to 300 L, with bacterial and cell versions plus great flexibility in configuration choices.

Professional

Validation of Techfors is a standard option, whereas compliance to GMP and other industrial standards can be assured as required.

Compact

Techfors is built on an "open ring frame" which takes up little space and even allows 150 L units to be easily moved from place to place.

Versatile

The range of options and choices is so great that almost any application can be accommodated and includes CIP (cleaning in place), magnetically coupled drives and swivel lids for easy access in confined spaces.

Modular Design

The Techfors allows almost any configuration, as for example:

- Bottom- or top-drive agitation system
- Classic stirrer shaft with mechanical seal or magnetic coupling
- Single or double mechanical seal
- Various types of stirrers and stirrer com-
- Manual or automatic sterilisation of the culture vessels
- Sterilisation of the air-inlet filter and peripherals
- Exhaust line with membrane or depth filter, heated with steam or electricity
- Heating via heat exchanger or direct injection of steam or hot water into the double jacket
- Combined harvest and sample-taking valve as well as separate valve systems by INFORS HT or by Nova Septic
- Membrane valves
- Connections for CIP
- A large range of accessories: spin filter, gas mix, etc.











Infors AG	Infors GmbH	Infors UK Ltd	Infors Sarl	Infors Benelux BV	Infors Canada	Infors Biotechnology
Headoffice, Switzerland						(China) Co., Ltd.
						Room 505C, Building 106
		The Courtyard Business	2, rue du Buisson	Fabriekstraat 38		Lize Zhongyuan
		Centre, Dovers Farm,	aux Fraises	NL-7005 AR Doetinchem	8350 rue Bombardier	Wangjing New Industrial Zone
Rittergasse 27	Dachauer Str. 6	Lonesome Lane, Reigate,	Bâtiment D13	P.O. Box 613	Anjou, Quebec	Chaoyang District, Beijing
CH-4103 Bottmingen	D-85254 Einsbach	Surrey, RH2 7QT, UK	FR-91300 Massy	NL-7000 AP Doetinchem	Canada H1J 1A6	100102 P. R. of China
T +41 (0)61 425 77 00	T +49 (0)8135 8333	T +44 (0)1737 22 31 00	T +33 (0)1 69 30 95 04	T +31 (0)314 364 450	T +1 514 352 50 95	T +86 (0)10 51 65 20 68
F +41 (0)61 425 77 01	F +49 (0)8135 8320	F +44 (0)1737 24 72 13	F +33 (0)1 69 30 95 05	F +31 (0)314 378 076	F +1 514 352 56 10	F +86 (0)10 64 39 05 85
headoffice@infors-ht.com	infors.de@infors-ht.com	infors.uk@infors-ht.com	infors.fr@infors-ht.com	infors.bnl@infors-ht.com	infors.ca@infors-ht.com	info@infors-ht.com.cn